

Six Channel Cordless Falls Monitor



PART NUMBER: 2700

FEATURES

- Modern ergonomic design
- Monitor up to 7 devices simultaneously - a single corded Pad / Mat and up to 6 wireless devices.
- Comes with batteries, durable ABS plastic enclosure and Caregiver Reset Key
- Alarm memory and auto re-arm
- Amber LED visual alarm indicator
- Caregiver secure reset option
- Battery operated or optional external plug pack
- Up to 6 months battery life (nominal)
- Easy to access battery compartment
- Can be connected to a nurse call system
- 3 selectable alarm tones & thumbwheel adjustable alarm volume
- "Pause" feature allows the caregiver to attend to a patient without setting the alarm off
- Alerts for lost nurse call connection, lost AC power supply, change pad (SafeTime) & low battery

SAFETY FEATURES

- Non resettable "Pad Signal Lost" warning
- Caregiver secure reset option
- All safety alerts sound even if alarm volume set to silent, including low battery, AC lost, pad signal lost, change pad (SafeTime) & nurse call lost
- Cable strain relief & breakaway cord technology prevents damage & helps reduce risk of entanglement when used with Premium Pads & Mats

TECHNICAL SPECIFICATIONS

- DC operation using either a 3 x AA batteries (included) or 6 V plug pack
- Dimensions: 120 x 77 x 40 mm
- Weight: 150 g without batteries
- Range: Up to 50 m* when used with cordless bed pads and chair pads
- Operating temperature: 5°C to 40°C
- Relative humidity: < 85%
- Operating frequency: 433 Mhz fully supervised
- Volt free normally open relay contact alarm output
- RJ10 type socket for optional wired bed/chair pad & floor mat connection
- Low battery detection
- Missing AC / DC power supply detection and indication

ACCESSORIES

All **Premium & Universal Pads & Mats**, All **Cordless Pads & Mats** (using a **2719** transmitter), **2406** Plug Pack 6 V, **2558** Wireless Call Button, **2595** Wristband Alarm Transmitter, **2588** Personal Emergency Transmitter (PET), **2599** Wireless SOS Call Button, **2636** Wheelchair / Chair Seatbelt, **2642** Caregiver reset key, **3500** ActiveCare Console, **3501** ActiveCare Pager, **3502** Cordless Handset, **3503** Wireless Door Contact, Multiple Nurse Call Cables

cura¹.com

*The range of wireless devices may vary and be dependent on your specific environmental conditions.
cura¹ products are designed to augment well designed care programs and practices, cura¹ does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision, and are at risk.